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DE | Deutsch

1. Produktbeschreibung

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EN | English

1. Product description

els cem is a light- and chemically-curing, radiopaque composite luting cement in a double-mix syringe. The color-stable material is offered in an universal transparent version. els cem is a polymer-based restorative material type 2, class 3 according to DIN EN ISO 4049 and has a radiopacity of 3.2 mm aluminium.

2. Composition
BisEMA, barium glass silanized, fumed silica, ytterbium trifluoride (Filler content: 61 % wt, 40 % by volume, particle size between 0,04 - 3,00 µm, median 0,7 µm) catalysts, inhibitors, pigments

3. Intended Use

Saremo cements and auxiliary bonding promoters are intended for the indirect reconstituted or correction of functionally compromised natural dentition (e.g., deficient teeth).

4. Indication

Permanent cementation of indirect composite, metal or porcelain/zirconia restorations (inlays/onlays, crowns, bridges, root posts).

5. Contra-indications

This cement is contraindicated for inlay-restorations, distooral prostheses, use on opened pulp or in case of pulpitis or known allergy to the above mentioned ingredients of the product.

6. Patient target group:

els cem is suitable for use in all patients without any age or gender restrictions.

7. User

els cem should only be used by a professionally trained dental practitioner.

8. Side effects

In individual cases, contact allergies have been described to products of a similar composition. In case of doubt, the allergy should be clarified and excluded by a specific test before the application of this product. In case of known hypersensitivity (allergies) to these ingredients of *els cem*, the application should be avoided. To avoid pulp reactions, it is recommended to expose the exposed dentine in the pulp area with suitable underfilling materials (preferably with calcium hydroxide).

Note: els cem does not contain BisGMA, TEGDMA or HEMA.

9. Interactions

Phenolic compounds (e.g., eugenol or preparations containing thymol), oxidative disinfectants (e.g., hydrogen per-oxide) or strong acids, in particular acidic, one-component adhesive systems or acid primers, can impair polymerisation of the product. Basic sandblasting abrasives inhibit the effect of the product SAREMCO *cmf prime*.

10. Processing stages

10.1. Preparing the restoration and tooth

Roughen (metal surfaces, zirconium oxide) respectively silanise (porcelain) inner surfaces of the restoration or workpiece to create a physical retention. Especially with zirconium oxide surfaces a good retentive pattern is decisive for the long-life cycle of the restoration. Blast the inner surfaces of zirconium oxide restorations with 110 µ aluminium oxide at 2 - 3 bar (roughen, enlarge surface) and then clean with oil-free air. Afterwards, the use of a zirconium oxide primer is recommended. Observe the manufacturers' instructions for use. Keep dry! Remove temporary restoration and any residues of the temporary restoration, clean and dry the core or cavity.

10.2. Enamel and dentine etching

Apply the etching gel and condition in correspondence to instruction for use of the product. SAREMCO recommends *cmf etch*. Sufficiently etched enamel appears matt and white.

10.3. Adhesion

Apply the adhesive and condition in correspondence to instruction for use of the product. For patients with a proven allergy to HEMA, TEGDMA or BisGMA SAREMCO recommends *cmf adhesive system* and *els duobond*.

Warning: To avoid interaction of this product with a strongly acidic medium, do not use any acidic one-component adhesives or cover acidic primers with a properly hardened adhesive film. Adhesives should always be applied in a very thin layer due to the fitting accuracy of the workpiece.

10.4. Cementation, fitting, curing

Remove the cap and place a new mixing tip. Dispense the product directly from the syringe onto the surfaces of the restoration to be bonded. Insert the restoration and fix it with slight pressure in the correct end position.

Warning: At a room temperature of 23°C / 73°F and a relative humidity of 50 %, the following times are available for the processing:

Processing time: > 1:30 min
Curing time: < 5:00 min

Lower temperatures extend the stated times, while higher temperatures or stronger exposure to light shorten them. The stated exposure time is based on light polymerization devices with an emission maximum ranging between 450 - 470 nm and a light power ≥ 700 mW/cm².

The necessary exposure time can vary depending on the light source and its instructions for use. The inserted working piece must not be moved after the end of the processing time.

Remove excess immediately by using a suitable instrument/aid. Light-cure the connections between tooth and rest-oration for 20 seconds (per segment). Please be aware that excess luting composite cement must be removed in time.

10.5. Removal of excessive cement

a) Wiping technique:
Remove excess cement immediately afterwards with a brush, floss or a scaler. In particular, pay attention on the timely removal of excess luting material in less accessible areas

(approximal, gingival margins, bridge pontic). If in doubt, please take an X-ray image.

b) Quarter technique (max. 6 bridge abutments) - light-curing of all excess material and subsequent removal:

The excess cement is to be light-activated using a curing light at a distance of max. 10 mm for 20 seconds per quarter side (mesiooral, distooral, mesiobukkal, distobuccal). This facilitates removal with a scaler. Again, pay attention to a timely removal of excess luting material in less accessible areas (approximal, gingival margins, bridge pontic). If in doubt, please take an X-ray image.

Carefully remove any excess cement in mucosal proximity. While removing excess cement, do not brush the cement out of the margins.

10.6. Finishing

Finish the margins of the restoration and proximal areas by using suitable instruments (finishing diamonds, flexible polishing discs, strips etc.). It is recommended to fluorideate the treated teeth.

11. Storage

Protect light-curing products from strong sources of light and heat! els cem was developed for use at room tempera-ture (20°C - 25°C / 68°F - 77°F). Store in the fridge between 4°C - 8°C / 39°F - 46°F. Do not freeze! Bring to room temper-ature before use.

12. Batch number and expiry date

The batch number should be specified to identify products in case of enquiries. Products should no longer be used once the expiry date has elapsed.

13. Precautionary measures

Close containers after each use with the right lid. Keep out of reach of children. For dental use only. Commercially available medical gloves do not provide protection against the sensitisation effect of methacrylates. If the product comes in contact with the glove, remove the glove and dispose of it, wash your hands with water and soap immediately and put on a new glove. In case of an allergic reaction, seek medical advice.

14. Emergency measures

In case of direct contact with the oral mucosa, rinse with water. In case of contact with the eyes, rinse thoroughly with water. Consult an eye specialist.

15. Hygiene

Use application instruments for one patient only. Dispose products away from patients to avoid contamination. Do not disinfect double syringes with oxidative disinfectants, rather with medical alcohol.

16. Warranty

The product was developed for use in dentistry and must be processed in accordance with the instructions for use. For further damages, namely that caused by non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. Our liability is restricted to the quality of our products. In the case of a product being of defective quality, only its value is replaced. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental AG.

17. Other notes for Europe

If the user and/or patient become aware of serious incidents connected with application of the product, they are to be reported to the manufacturer and the responsible authorities of the state in which the user and/or patient resides. The Summary of Safety and Clinical Performance of SAREMCO Products can be found in the European database on medical devices (EUDAMED-<https://ec.europa.eu/tools/eu-damed/>).

18. Production / distribution
SAREMCO Dental AG
Gewerbestrasse 4
CH-9445 Rebstein / Switzerland
Tel: +41 (0) 71 775 80 90
Fax: +41 (0) 71 775 80 99
info@saremo.ch
www.saremo.ch

Edited 07-2024 | D600137

Class IIa medical devices
Class III medical devices (Canada)

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