# SAREMCO

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# SAREMCO

## Instruction for use james-2

**1. Product description**  
Light-curing, solvent-free 1-component restorative adhesive in the total-etch process for enamel and dentine.

**2. Composition**  
hydroxyethylmethacrylate, urethanemethacrylate, polyalkenatoate methacrylated, hydroxypropylmethacrylate, glycerinedimethacrylate, catalysts, inhibitors

**3. Intended Use**  
Saremco dental bonding agents promote adhesion of direct or indirect prosthesis for reconstruction or correction of functionally compromised natural dentition (e.g., deficient teeth).

**4. Indication**  
1. Direct restorations of all cavity classes on anterior and posterior teeth with light-curing composites.  
2. Direct restorations and core build-ups with dual-curing or self-curing composites  
3. Adhesive bonding of indirect restorations (crowns, bridges, inlays, onlays, veneers) with light-curing, dual-curing or self-curing composite cements

**5. Contra-indication**  
Opened pulp, pulpitis, known allergy to methacrylates.

**6. Patient target group**  
james-2 is suitable for use in all patients without any age or gender restrictions.

**7. User**  
james-2 should only be used by a professionally trained dental practitioner.

**8. Side effects**  
In individual cases, contact allergies have been described to products of a similar composition. To avoid pulp reactions, it is recommended to cover the exposed dentine in the pulp area with suitable underfilling materials (preferably with calcium hydroxide).

**9. Interactions**  
Avoid underfilling materials which may hinder polymerisation due to their ingredients. All phenolic compounds, such as zinc oxide eugenol or preparations containing thymol, belong to this category.

### 10. Processing stages

**10.1. Drying**  
Due to the following application of the filling material, a drainage of the working field is required. Avoid any contamination of the cavity with blood or saliva. A rubber dam is recommended.

**10.2. Cavity preparation**  
Prepare the cavity as usual. Undercuts and bevelled margins are recommended to improve adhesion and the margin fit. Clean and dry the cavity. Avoid overdrilling.

**10.3. Etching of enamel and dentine**  
Apply the etching gel and condition in correspondence to instruction for use of the product. SAREMCO recommends *microd* and *cmf etch*. Sufficiently etched enamel appears matt and white.

**Note:** Possibly use the total-etch-technique in cavities limited by dentine. In this case, the dentine should not be etched for more than 15 seconds and only dried briefly after rinsing.

**10.4. Application**  
Apply *james-2* to enamel/dentine in a not too thin layer with a brush and massage into the tooth substance for 20 seconds.

**Note:** *james-2* is light-curing, so avoid too intense ambient light. Dim the surgical light during the application.

**8.5. Drying**  
Dry the adhesive layer with an air blower for at least 5 seconds.

**8.6. Polymerizing**  
Polymerize the bonding layer with blue light (halogen or LED lamp) for 30 seconds. The specified exposure time relates to halogen or LED light-curing devices with a minimum light intensity of 500 mW/cm<sup>2</sup> and a wavelength of 400 - 500 nm.

**Note:** Do not use *james-2* for opaque restorations or deep cavities where complete light-curing is not guaranteed.

Application and processing of the composite filling material according to the instructions for use of the product used. SAREMCO recommends *mhc composite*, *els composite* and *apt composite*.

**11. Storage**  
Close *james-2* bottle tightly after use. Store in an upright position to ensure the backflow of the liquid. Do not expose light-curing products to sunlight or operating light. *james-2* was developed for use at room temperature (20°C - 25°C / 68°F - 77°F). Store at temperatures between 4°C - 28°C / 39°F - 82°F. If room temperature exceeds 28°C / 82°F storing in the refrigerator is recommended. Do not freeze! Ambient temperatures above 28°C / 82°F can reduce the shelf-life of the product.

**12. Batch number and expiry date**  
The batch number should be specified to identify products in the case of enquiries. Products should no longer be used once the expiry date has elapsed.

**13. Precautionary measures**  
Close containers after each use with the right lid. Keep out of reach of children. For dental use only. Commercially available medical gloves do not provide protection against the sensitisation effect of methacrylates. If the product comes into contact with the glove, remove the glove and dispose of it, wash your hands with water and soap immediately and put on a new glove. In case of an allergic reaction, seek medical advice.

**14. Emergency measures**  
In case of direct contact with the oral mucosa, rinse with water. In case of contact with the eyes, rinse thoroughly with water. Consult an eye specialist.

**15. Hygiene**  
Use application instruments for one patient only. Dose products away from patients to avoid contamination.

**16. Warranty**  
The product was developed for use in dentistry and must be processed in accordance with the instructions for use. For further damages, namely that caused by non-compliance with the instructions for use or other improper handling or inappropriate use of a product, any liability is rejected. Our liability is restricted to the quality of our products. In the case of a product being of defective quality, only its price is replaced. It is the responsibility of the user to check, before using the products, whether they are suitable for the intended purpose. He expressly assumes all risks associated with using the product and is solely responsible for any resulting damages. Safety data sheets and technical data sheets are available on the website of SAREMCO Dental AG.

**17. Other notes for Europe**  
If the user and/or patient becomes aware of serious incidents in connection with the use of the product, they are to be reported to the manufacturer and the responsible authorities of the state in which the user and/or patient resides. The Summary of Safety and Clinical Performance of SAREMCO Products can be found in the European database on medical devices (EUDAMED-<https://ec.europa.eu/tools/eudamed/>).

**18. Production / distribution**  
SAREMCO Dental AG  
Gewerbestrasse 4  
CH-9445 Rebstein / Switzerland  
Tel: +41 (0) 71 775 80 99  
Fax: +41 (0) 71 775 80 99  
info@saremco.ch  
www.saremco.ch

Edited 07-2024 | D600170

Class IIa medical devices  
Class III medical devices (Canada)

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